

Percent Moisture

Method SOP642

Lab Name: Paragon Analytics

Date Extracted: 06/03/2004

Date Analyzed: 06/03/2004

Analyst: Crystal Halverson

Validated By: ckh

Validation Date: 06/03/2004

Validation Time: 9:55:09 AM

Run ID	Prep Batch ID	QC Batch ID	Lab ID	QC Type	Dish Wt	Wet Wt	Dry Wt	Dry Wt-Dish Wt	Percent Moisture	Percent Solids	RPD
EX040602-4A	EX040602-4	EX040602-4-1	0405215-1	SMP	1.35	11.19	12.13	10.78	3.7	96.3	
EX040602-4A	EX040602-4	EX040602-4-1	0405215-2	SMP	1.35	10.32	11.2	9.85	4.6	95.4	
EX040602-4A	EX040602-4	EX040602-4-1	0405215-3	DUP	1.36	10.3	11.06	9.70	5.8	94.2	
EX040602-4A	EX040602-4	EX040602-4-1	0405215-3	SMP	1.35	11.94	12.53	11.18	6.4	93.6	

QC Types

CAR	Carrier reference sample	DUP	Laboratory Duplicate
LCS	Laboratory Control Sample	LCSD	Laboratory Control Sample Duplicate
MB	Method Blank	MS	Laboratory Matrix Spike
MSD	Laboratory Matrix Spike Duplicate	REP	Sample replicate
SMP	Field Sample	SYS	Sample Yield Spike

Comments:

DUP = Sample Duplicate

Wet Wt = Sample Wet Wt - Dish Wt

Dry Wt = Sample Dry Wt + Dish Wt

Dry Wt - Dish Wt = Sample Dry Wt - Dish Wt

All weight values shown above are expressed in grams.